

Custom trays

Selecting your impression tray size.


The following template has been provided to assist you in selecting the most appropriate impression tray applicable to your mouth size.

Selecting the most appropriate impression tray assists in manufacturing the best possible professional custom whitening tray.

- ① Print out this template.
- ② With your index finger and your thumb measure approximately midway across your upper mouth.
- ③ With your fingers open at the width of your upper mouth, place your fingers against each impression tray printed on the template to determine the most appropriate size for your mouth.
- ④ View and print the instructions here on how to make the impression for custom fitted tooth whitening trays.

Impression tray templates

Large impression Tray (8cm)



Medium impression Tray (7cm)



Small impression Tray (6cm)

